RadiForce® MX240W

Color LCD Monitor

English

Deutsch

Français

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Instructions for Use: ■ Setup Manual

Important: Please read this Setup Manual and the User's Manual stored on the CD-ROM carefully to familiarize yourself with safe and effective usage. Please retain this manual for future reference.

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
Color LCD Monitor

SAFETY SYMBOLS

This manual uses the safety symbols below. They denote critical information. Please read them carefully.

 WARNING Failure to abide by the information in a WARNING may result in serious injury and can be life threatening.	 CAUTION Failure to abide by the information in a CAUTION may result in moderate injury and/or property or product damage.
 Indicates a prohibited action.	
 Indicates a mandatory action that must be followed. For example, the symbol  illustrates the notification of general prohibition such as "Grounding the unit".	

PRECAUTIONS

WARNING

If the unit begins to emit smoke, smells like something is burning, or makes strange noises, disconnect all power connections immediately and contact your dealer for advice.

Attempting to use a malfunctioning unit may result in fire, electric shock, or equipment damage.

Use the enclosed power cord and connect to the standard power outlet of your country.

Be sure to remain within the rated voltage of the power cord. Not doing so may result in fire or electric shock.

"Power Supply: 100-120/200-240 Vac 50/60Hz"

Keep small objects or liquids away from the unit.

Small objects accidentally falling through the ventilation slots into the cabinet or spillage into the cabinet may result in fire, electric shock, or equipment damage. If an object or liquid falls/spills into the cabinet, unplug the unit immediately. Have the unit checked by a qualified service engineer before using it again.

Use the unit in an appropriate location.

Not doing so may result in fire, electric shock, or equipment damage.

- Do not place outdoors.
- Do not place in the transportation system (ship, aircraft, trains, automobiles, etc.).
- Do not place in a dusty or humid environment.
- Do not place in a location where water is splashed on the screen (bathroom, kitchen, etc.)
- Do not place in a location where the steam comes directly on the screen.
- Do not place near heat generating devices or a humidifier.
- Do not place in a location where the product is subject to direct sunlight.
- Do not place in an inflammable gas environment.

The equipment must be connected to a grounded main outlet.

Not doing so may result in fire or electric shock.

Fluorescent backlight lamps contain mercury (the products that have LED backlight lamps contain no mercury), dispose according to local, state or federal laws.

Exposure to elemental mercury can result in effects on the nervous system, including tremor, memory loss, and headache.

CAUTION

Do not block the ventilation slots on the cabinet.

- Do not place any objects on the ventilation slots.
 - Do not install the unit in a closed space.
 - Do not use the unit laid down or upside down.
- Blocking the ventilation slots prevents proper airflow and may result in fire, electric shock, or equipment damage.

Use an easily accessible power outlet.

This will ensure that you can disconnect the power quickly in case of a problem.


Package Contents

Check that all the following items are included in the packaging box. If any items are missing or damaged, contact your local dealer.


NOTE

Please keep the packaging box and materials for future movement or transport of the monitor.


- Monitor
- Power cord
- Digital signal cable : FD-C39

- Analog signal cable : FD-C16


- EIZO USB cable : MD-C93

- EIZO LCD Utility Disk (CD-ROM)
 - RadiCS LE (for Windows)
 - ScreenManager Pro for Medical (for Windows)
 - User's Manual
- Setup Manual (this manual)
- Limited warranty
- Recycling Information

Controls and Functions

1	Sensor	Detects ambient brightness. BrightRegulator function. Refer to the User's Manual on the CD-ROM.
2	Adjustment Lock button	This function locks the buttons to retain the status adjusted or set once.
3	Input Signal Selection button	Switches input signals for display when two PCs are connected to the monitor.
4	Mode button	Displays the CAL Switch menu (page 5).
5	Auto button	Performs the function to adjust the screen automatically. (analog input only)
6	Enter button	Displays the Adjustment menu, determines an item on the menu screen, and saves values adjusted.
7	Control buttons (Left, Down, Up, Right)	Chooses an adjustment item or increases/decreases adjusted values for advanced adjustments using the Adjustment menu.
8	Power button	Turns the power on or off.
9	Power indicator	Indicates monitor's operation status. Blue: Operating Orange: Power saving Off: Power off
10	Security lock slot	Complies with Kensington's MicroSaver security system.
11	Power connector	Connects the power cord.
12	Input signal connectors	DVI-I Connector x 2
13	USB port (Up)	Connects the USB cable to use the software that needs USB connection, or to use USB Hub function.
14	USB port (Down)	Connects a peripheral USB device.
15	Stand	Used to adjust the height and angle of the monitor screen.
16	Cable holder	Covers the monitor cables.

* ScreenManager® is an EIZO's nickname of the Adjustment menu. (For how to use ScreenManager, refer to the User's Manual on the CD-ROM.)


Connecting Cables

Note

- Check that the monitor and the PC are powered off.
- When replacing the current monitor with a MX240W monitor, be sure to change the PC settings for resolution and vertical frequency to those which are available for the MX240W monitor referring to the resolution table (back of cover page) before connecting the PC.

Tips

- When connecting two PCs to the monitor, refer to the User's Manual on the CD-ROM.

- 1 **Connect the monitor to the PC with a signal cable that matches the connectors.**
After connecting the cable connectors, tighten the screws of the connectors to secure the coupling.

Digital Input

Connectors on the monitor	cable	Connectors on the PC
	Signal cable FD-C39 	


Analog Input

Connectors on the monitor	cable	Connectors on the PC
	Signal cable FD-C16 	

- 2 **Connect the USB cable when using RadiCS LE or ScreenManager Pro for Medical.**

Connecting to the Monitor

Connecting to a Computer

- 3 **Plug the power cord into a power outlet and the power connector on the monitor.**

- 4 **Press to turn on the monitor.**

The monitor's power indicator lights up blue.

- 5 **Turn on the PC.**

The screen image appears.

If an image does not appear, refer to "No-Picture Problem" (page.6) for additional advice.

- 6 **When using analog input signals, perform the Auto Adjustment function (Refer to "Performing Auto Adjustment function" below).**


The Auto Adjustment function is not necessary when digital signals are input because images are displayed correctly based on the preset data of the monitor.

Note

- Turn off the monitor and PC after using them.
 - For the maximum power saving, it is recommended that the Power button be turned off. Unplugging the power cord completely shuts off power supply to the monitor.
-

Adjusting the Screen Height and Angle

Hold left and right edge of the monitor with both hands, and adjust the screen height, tilt and swivel of the screen to the best condition for working.

Note

- Be sure that the cables are correctly connected.
-

Installing RadiCS LE

The quality control software "RadiCS LE" (for Windows) is included with MX series monitors. To calibrate the monitor and manage the calibration history, install RadiCS LE with reference to the RadiCS LE User's Manual on the CD-ROM.

Note

- Connect the monitor to the PC with a USB cable when using RadiCS LE.
-

Tips

- To perform advanced quality control that conforms to medical standards, use the optional monitor quality control software kit "RadiCS UX1".
-

Performing Auto Adjustment Function (Analog input only)

When analog signals are input, the Auto Adjustment function enables the automatic adjustment of clock, phase, screen position, and resolution. For details of the Auto Adjustment function, refer to the User's Manual on the CD-ROM.

1 Press **ⓐ**.

A message window appears.

2 Press **ⓐ** again while the message is displayed.

The Auto Adjustment function is activated and clock, phase, and screen position are adjusted automatically.

Message window

Your setting will
be lost if you
press again now.

Note

- The Auto Adjustment function works correctly when an image is fully displayed over the display area of a Macintosh or Windows PC screen. It does not work properly when an image is displayed only on a part of the screen (DOS prompt window, for example) or when a black background (wallpaper, etc.) is in use.
 - The Auto Adjustment function may not work properly with some graphics boards.
-

Selecting Display Mode

CAL Switch function allows you to select the best display mode easily according to monitor's application.

For details of CAL Switch function, refer to the User's Manual on the CD-ROM.

CAL Switch Mode

DICOM-CL	Available for setting tailored to X-ray film color (clear base).
Text	Suitable for displaying texts for word processing or spreadsheets.
Custom	Available for making desired setting.
CAL	Used for monitor calibration.

1 Press .

The CAL Switch menu appears.


2 Press  again while the CAL Switch menu is displayed.

Each time the button is pressed, the mode on the screen switches. (Refer to the CAL Switch Mode table.)

3 Press  at the desired mode.

The selected mode is set.

CAL Switch menu

The basic settings and adjustments are completed. For advanced settings/adjustments, refer to the User's Manual on the CD-ROM.

No-Picture Problem



If no picture is displayed on the monitor even after the following remedial action is taken, contact your local dealer.

1. Check the Power indicator.

Symptom	Status	Possible cause and remedy
No picture	Power indicator does not light up.	<ul style="list-style-type: none"> • Check whether the power cord is connected correctly. If the problem persists, turn off the monitor, and then turn it on again a few minutes later. • Press .
	Power indicator lights blue.	<ul style="list-style-type: none"> • Set each RGB adjusting value in <Gain> to higher level.
	Power indicator lights orange.	<ul style="list-style-type: none"> • Switch the input signal with . • Press a key on the keyboard or click the mouse. • Turn on the PC.

2. Check the error message that remains on the screen.

These messages appear when the input signal is incorrect even if the monitor is functioning.

Symptom	Status	Possible cause and remedy
 <p>Signal Check</p> <p>Signal 2</p> <p>fH: 0.0kHz</p> <p>fU: 0.0Hz</p>	The input signal is not received.	<ul style="list-style-type: none"> • Check whether the PC is turned on. • Check whether the signal cable is properly connected to the PC. • Switch the input signal with .
 <p>Signal Error</p> <p>Signal 2</p> <p>fD: 165.0MHz</p> <p>fH: 75.0kHz</p> <p>fU: 60.0Hz</p>	The signal frequency is outside the specification. Incorrect signal frequency is shown in red.	<ul style="list-style-type: none"> • Reboot the PC. • Change the mode to an appropriate mode using the graphics board's utility software. Refer to the manual of the graphics board for details.

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No part of this manual may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, or otherwise, without the prior written permission of EIZO Corporation.

EIZO Corporation is under no obligation to hold any submitted material or information confidential unless prior arrangements are made pursuant to EIZO Corporation's receipt of said information. Although every effort has been made to ensure that this manual provides up-to-date information, please note that EIZO monitor specifications are subject to change without notice.

Apple, Mac, Macintosh, iMac, eMac, Mac OS, MacBook, PowerBook, ColorSync, QuickTime, and iBook are registered trademarks of Apple Inc.

Windows, Windows Vista, Windows Media, SQL Server, and Xbox 360 are registered trademarks of Microsoft Corporation in the United States and other countries.

VESA is a registered trademark or a trademark of the Video Electronics Standards Association in the United States and other countries.

EIZO, the EIZO Logo, ColorEdge, DuraVision, FlexScan, FORIS, RadiForce, RadiCS, RadiNET, Raptor, and ScreenManager are registered trademarks of EIZO Corporation in Japan and other countries.

All other company and product names are trademarks or registered trademarks of their respective owners.

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<http://www.eizo.com>

Copyright © 2009-2014 EIZO Corporation All rights reserved.

4th Edition-January, 2014 Printed in Japan.

00NOL505D1
(U.M-SURMX240W-4)

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EIZO Europe GmbH